

# Wireless medical keyboard Easy to clean and disinfect



Purekeys medical keyboards combine a completely flat top with the tactile feedback of a normal keyboard. The flat top surface is made of a high grade silicone rubber and can be cleaned with any

hospital grade disinfectant. With cleaning time less than 5 seconds, best practice to clean between patients can be easily achieved.


Flat surface: easy to clean without the need for immersion, no raised surfaces where bacteria can hide

Real scissor type keys with excellent tactile key response

Works with any Microsoft Windows-based PC

Silicone based cover, non allergenic, can be cleaned with hospital grade antibacterial fluids

Dedicated Clean Key to quickly disable input during cleaning procedures



L/W: 457x148 mm  
Weight incl. 4 x AAA batteries: 950 grams  
IP66, can be cleaned with liquid disinfectants  
Color: White with black and blue print  
Interface: RF 2.4 GHz, easy installation, dongle already paired  
Battery life: between 2 and 4 years

Approvals: CE, FCC  
Warranty: 2 Year limited Warranty  
Available in the following layouts: US-Int, Scandinavian, UK, Arabic, German, French, Spanish, Italian, Hebrew, Portuguese, Belgium, Swiss.  
Other layouts available upon request.  
part# 500XXX00, XXX=country code for layout

Also available from Purekeys:  
Disinfectable wireless mouse  
with touch scroll.  
part# 40070



 **PUREKEYS**  
Medical Keyboards

Purekeys BV  
Rouaanstraat 23 C  
9723 CC Groningen  
The Netherlands

Tel. +31 507 210 810  
Fax. +31 507 210 811  
E-mail [info@purekeys.com](mailto:info@purekeys.com)  
[www.purekeys.com](http://www.purekeys.com)