

**ORIGINAL AIR-FLOW®  
THERAPY → SUPRAGINGIVAL**

**CLINICAL APPLICATIONS**

**AIR-FLOW®  
HANDPIECES**



**AIR-FLOW®** > Sodium bicarbonate  
**POWDER** > Average grain size ~65µm  
**CLASSIC** > Available in 6 different tastes

300 g bottle

- REMOVAL OF BIOFILM AND HEAVY DISCOLORATIONS
- SURFACE PREPARATION BEFORE DENTAL TREATMENTS



**AIR-FLOW®** > Sodium bicarbonate  
**POWDER** > Average grain size ~40µm  
**CLASSIC**  
**COMFORT**

250 g bottle

- REMOVAL OF BIOFILM AND HEAVY DISCOLORATIONS  
WITH **OPTIMIZED PATIENT COMFORT**
- SURFACE PREPARATION BEFORE DENTAL TREATMENTS



**AIR-FLOW®** > Glycine-based  
**POWDER** > Average grain size ~65µm  
**SOFT**

200 g bottle

- REMOVAL OF BIOFILM AND MEDIUM/LIGHT DISCOLORATIONS
- TOOTH POLISHING



AIR-FLOW®

**→ SUB- AND SUPRAGINGIVAL**



**AIR-FLOW®** > Erythritol-based  
**POWDER** > Average grain size ~14µm  
**PLUS** > Contains chlorhexidine (0,3%)

100 g bottle

- REMOVAL OF BIOFILM AND MEDIUM/LIGHT DISCOLORATIONS  
WITH **OPTIMIZED PATIENT COMFORT**
- TOOTH POLISHING
- REMOVAL OF BIOFILM FROM POCKETS **DOWN TO 4mm**  
ON NATURAL TEETH AND IMPLANTS
- REMOVAL OF BIOFILM FROM DEEP PERIODONTAL POCKETS,  
**FROM 4 TO 10mm** ON NATURAL TEETH AND IMPLANTS



AIR-FLOW®  
AIR-FLOW® PLUS

**→ SUBGINGIVAL**



**AIR-FLOW®** > Glycine-based  
**POWDER** > Average grain size ~25µm  
**PERIO**

120 g bottle

- REMOVAL OF BIOFILM FROM POCKETS **DOWN TO 4mm**  
ON NATURAL TEETH AND IMPLANTS
- REMOVAL OF BIOFILM FROM DEEP PERIODONTAL POCKETS,  
**FROM 4 TO 10mm** ON NATURAL TEETH AND IMPLANTS



PERIO-FLOW®  
AIR-FLOW® PERIO